

STEP 1 - GETTING STARTED: DETAILED PROGRESS CHECKLIST



LIST OF REVISIONS

REVISION No.	DATE	REVISED BY:	Page/s	DESCRIPTION
1	TBD	GDSP	4	Draft Version





DETAILED PROGRESS CHECKLIST

The purpose of this document is to provide detailed activities that HIE partner needs to perform to complete the onboarding process. These are list of all the activities that should be completed during the implementation phase, interface testing phase and for production go-live.

Phase	Item	Status
Discontinu	Initiate project, identify resources and conduct project kick-off meeting	
Planning	Schedule weekly status meeting with GDSP	
	Identify the scope of implementation	
	Complete the understanding of the HL7 message structure and standards used by GDSP	
	Complete the understanding of different panels used in HL7 NBS result messages	
Design	Outline the process to match the HL7 result with the medical record in the system using the identifiers from TRF forms	
	Outline the process to download result message zip from FTP folder	
	Identify the fields that will be used to handle twin birth results	
	Complete development to parse all segments (like MSH, PID, NK1, ORC, OBR and OBX) from HL7 messages	
Development	Complete development to extract information from narrative summary OBX from OBR 2	
20.0.0	Complete development to zip and unzip the result and acknowledgment messages	
	Complete development of error queue for unmatched results	



STEP 1 - GETTING STARTED: DETAILED PROGRESS CHECKLIST



Development	Complete development to generate positive (AA-type) Acknowledgement	
	Complete development to generate negative acknowledgement (AE and AR type)	
	Complete development to parse Additional Results as per HIE Message Generator Tool	
	Complete development to generate hI7 Acknowledgment files and zip with recommended file name format	
Testing	Complete HL7 result parsing for all types of messages generated by GDSP using HIE Message Generator Tool	
	Complete validating additional results using HIE HL7 generator tool	
	Complete validating positive (AA type) acknowledgement using HIE Validator tool	
	Complete validating negative (AE & AR type) acknowledgement using HIE Validator tool	
Review	Review the screening result report with your organization Medical Director	
	Send the screening result report to GDSP	
	Confirm approval from Medical Director on the screening result report format	
	Inform stakeholders for Interface Testing	
Interface	Identify internal resources for Interface testing	
Testing	Identify GDSP resources for Interface testing	
	Setup infrastructure for Interface testing	



STEP 1 - GETTING STARTED: DETAILED PROGRESS CHECKLIST



	Complete Interface testing with GDSP	
	Inform stakeholders for Production go-live	
	Identify internal resources for go-live	
Go-Live	Identify GDSP resources for go-live	
	Setup infrastructure for go-live	
	Send signed Consent forms to GDSP	